


# VG 55 Ventilator

## High-Flow Therapy

- An air/oxygen blender, allowing for 21% to 100% FiO<sub>2</sub>, generates up to 60L/min flow.



## Complete Ventilation Modes

- Volume Guaranteed Pressure Support Ventilation (VGPS) assures tidal volume with auto adjusting pressure support, decreasing the workload for clinicians while ensuring patient safety.
- CPAP with Crelief option offers two levels of flow-based expiratory pressure relief, which improves patient comfort and treatment acceptance.



Corporate Headquarter  
HEYER Medical AG

Carl-Heyer-Str. 1/3, 56130 Bad Ems, Germany  
Tel: +49 2603 791 0  
Fax: +49 2603 791 200  
info@heyarmed.com  
www.heyarmed.com

Manufacturer  
HEYER Medical Co., Ltd.

No.10, Chaobai Street, Yanjiao Development Zone  
065201, Sanhe City, Hebei, China  
Tel: +86 10 5841 1198  
Fax: +86 10 6371 8989  
marketing@heyarmed.com

The \* marked accessories/ functions are optional.

VG55-2024.01

CE 0123

# VG 55 Ventilator



# VG 55

## Ventilator

### Start Your Natural Breath with VG55

VG55 is specially designed for Non-Invasive Ventilation (NIV). Research has shown that early application of NIV can effectively decrease intubation rates and hospitalization length. The application of NIV can also prevent ventilator associated pneumonia (VAP), as well as promote respiratory care. VG55 adopts advanced turbine technology; provides NIV with ICU-grade monitoring.



15 Seconds

#### Start Your Natural Breath within 15 Seconds

- Auto-calibration is performed within 15 seconds, safely guaranteeing the ventilator's performance.
- Quick setup with guide-way system and intuitive interface.
- Smart trigger combines pressure and flow triggers, highly improving the patient-machine synchrony, allowing the patient to breathe in a more natural manner.

5 Minutes

#### Become a Respiratory Therapist in 5 Minutes

- With 15.6" intelligent touchscreen, the VG55 displays complete ventilation monitoring in a single screen.
- Simplifies the patient management with user-friendly interface and guide-way information.
- 3000 pieces event/alarm logs provide clear information to ensure continuity of patient care.


in Seconds

#### Get What You Need in Seconds


- Two minutes of emergency oxygen therapy provides 100% oxygen for whenever and wherever needed.
- Extended SpO<sub>2</sub>, PR module fulfill clinical requirements.\*
- High-Flow Therapy up to 60L/min.




Standby mode for patient-clinician interaction



Various information interfaces: USB, VGA, Nurse Call



Smart trigger and cycling ensures patient-machine synchrony



Alarm reset helps maximize patient safety while reducing excessive alarms